




RCS MEDIA



Product List (for use with RCS, RCS+, RCS Isolator, RCS HighFlow)

Product No	Product	Description
941 105	Agar Strips TC	Tryptic Soy Agar for Total Count
941 115	Agar Strips TC gamma irradiated	Tryptic Soy Agar for Total Count double wrapped
941 125	Agar Strips TCI gamma irradiated	Tryptic Soy Agar with H2O2 neutralisers, double wrapped
941 135	Agar Strips TSM	Modified Tryptic Soy Agar for total count
941 195	Agar Strips YM	Rose Bengal Agar for Yeasts and Moulds
941 205	Agar Strips SDX	Sabouraud Dextrose Agar for Yeasts and Moulds
941 215	Agar Strips SDX gamma irradiated	SDX Agar for yeast and moulds double wrapped
941 225 025	Agar Strips DG18	Dichloran Glycerol Agar for yeasts and moulds
941 405 025	Agar Strips S	Manitol Salt Agar for Staphylococci
941 505 025	Agar Strips C	MacConkey Agar for Coliform bacteria
941 705	Agar Strips Penase	For total count in air containing penicillin
941 715	Agar Strips LAC	For total count in air containing beta -lactams
941 605	Blank Strip Kit	Empty sterile agar strips
941 900	Plastic Cover Slide	To seal the agar strips